


# Comparing Augmentative and Alternative Communication (AAC) Systems

	Speak for Yourself				
					
Specifications					
Hardware/operating system requirements	iPad running iOS 8.1 and above				
Number of pre-programmed words out of the box	4,994				
Vocabulary capacity (maximum number of buttons)	14,000				
AAC language systems available	Speak for Yourself				
Navigational demands: Maximum number of touches to access vocabulary	2				
Number of button size/overlay configurations	1				
Symbols used and size of symbol library	Smarty Symbols® 11,000				
Voice options available	10 voices included Cereproc/Apple Acapela voices available for \$24.99 USD in app purchase				
Access methods	Direct selection				

# Comparing Augmentative and Alternative Communication (AAC) Systems



Keyguard available	Yes				
Language(s) available	English				
Environmental controls	No				
Ability to hide/mask buttons	Yes				
Price	\$199.99 USD				
<b>Vocabulary Organization</b>					
Core vocabulary	Yes				
Motor learning/memory	Yes. Every button has a unique motor plan.				
Categories	Yes, organized under each core word				
QWERTY keyboard	Yes with ability to install third party keyboards				
Single Meaning Symbols	Yes				
Visual Scenes	No				
Activity (Situation) Based	No				
<b>Features</b>					
Ability to search for words within vocabulary	Multi sensory search feature				
Backup/sharing of student's vocabulary data file	iTunes, Dropbox, email				


# Comparing Augmentative and Alternative Communication (AAC) Systems



Ability to share contents of message window	text, email, copy/paste, Twitter				
Data logging/tracking	History feature with summary within app				
Unique features	Babble Hold That Thought Whisper Mode				
<b>Editing/ Programming</b>					
Computer Programming Option	No				
Minimum number of steps (touches) to add a word to the system	6				
Ability to change: Button background color Text color Pronunciation Symbol	yes no yes yes				
Number of steps (touches) to mask/hide buttons	3 (with ability to hide many buttons at once)				
<b>Support</b>					
Website	<a href="http://www.speakforyourself.org">www.speakforyourself.org</a>				
Facebook	Speak for Yourself AAC				
Facebook User's Group	Speak for Yourself User's Group				

# Comparing Augmentative and Alternative Communication (AAC) Systems



					
Twitter	@Speak4AAC				
email	<a href="mailto:speakforyourselfaac@yahoo.com">speakforyourselfaac@yahoo.com</a>				
phone	215-605-0508				
YouTube	SpeakforYourselfapp				
Free codes provided to SLPs/evaluators	Yes				

Created by Speak for Yourself LLC 2016 Other AAC developers and manufactures are invited to add their information. Please contact Heidi LoStracco, MS, CCC-SLP for an editable chart.

V 1.0 10/12/16